


*The New Movement In Orthodontics™*

# NeoLucent Plus Ceramic Bracket System



# NeoLucent Plus Ceramic Bracket System

The *NeoLucent® Plus™* Low-Profile Ceramic Bracket System delivers superb aesthetics by blending with the natural tooth tone for a virtually invisible appearance. With improved performance due to a *Crunch COAT™* Base, more anatomically contoured base radii, and improved materials for a smoother finish and color consistency, *NeoLucent Plus* Ceramic Brackets provide the definitive performance you demand.

*Henry Schein® Orthodontics™* is dedicated to providing world class products and solutions to meet your every orthodontic need. As a full-service partner, we offer the tools to help expand your practice and diversify treatment options. So when your patients desire a cosmetic alternative, we have the right choice for you, and them.

## WHAT'S IN THE PLUS?

*Crunch Coat* Base

**Base surface is proprietary to Henry Schein Orthodontics**

**Bond Strength Comparable to Mesh Pad**

Fine alumina particles bonded to the bracket base provide a large surface area for adhesive bonding

**Consistent Debonding**

Because debonding occurs near the *Crunch COAT* adhesive interface, there is minimal distress to the enamel surface

## PATIENT SATISFACTION

**Invisible Appearance**

Crafted to the highest aesthetic standards

**Translucent, natural appearance** from

99.9% pure Polycrystalline Alumina

**Stain resistant**

throughout treatment



## COMFORT

**Designed with your patient's comfort in mind**

**Low profile**

**Smooth, polished surface**

**Rounded Corners**



SEM Picture

## COMPATIBLE IN/OUT WITH MAESTRO LOW-PROFILE BRACKET SYSTEM AND NEOCRISTAL SAPPHIRE BRACKET SYSTEM

For aesthetic anterior and  
metal posterior bracketing

### CLINICAL PERFORMANCE

**Accurate placement**  
achieved with an easily  
removable color-coded  
ID mark


**Strong bonding**  
and predictable  
debonding with a  
*Crunch COAT Base*

**Superior Sliding  
Mechanics** via  
a precision-cut  
archwire slot


**Non-allergenic** for  
sensitive patients as  
brackets do not contain  
nickel, chromium or  
other metals

**Easy Ligation** with ample  
tie wing clearance


**Routine Protocol** with  
no special instruments  
or adhesives required



**Maestro®** Bracket



**NeoLucent Plus** Bracket



**NeoCrystal™** Bracket

### SINGLE HOOK

**More Comfort for the Patient,**  
less bulky than universal hooks



The mission of *Henry Schein Orthodontics* is to provide state-of-the-art orthodontic products and innovative clinical solutions that enable our customers to offer exceptional patient care, while expanding the scope and profitability of their practices.



*The New Movement In Orthodontics™*

For more information on our products and educational offerings, please contact us:  
**HenryScheinOrtho.com** • **Customer Service:** 888.851.0533 • **In the U.S.:** 800.547.2000  
**Outside the U.S.:** +(1) 760 448 8600 • **Canada:** CERUM 800.661.9567

**To email an order:** USASales@HenryScheinOrtho.com  
**To email an order outside the U.S.:** IntlOrders@OrthoOrganizers.com  
Ortho Organizers, Inc. 1822 Aston Ave. Carlsbad, CA 92008-7306 USA

