



A collage of various dental and medical products. In the foreground, there are stacks of white disposable toothbrushes in a box labeled 'TOPBRUSH', several white cold cushions in a box labeled 'AKZENTA TOP FROST', a box of orange nitrile gloves labeled 'NON LATEX MASK', and a box of blue and white paper cups. There are also stacks of blue and white paper napkins, a roll of pink paper, and a roll of green paper. The background is a scenic view of a snowy mountain range under a blue sky with clouds. The text 'DENTAL D' and 'PRODUCT' is visible in the bottom right corner.

TABLE OF CONTENTS

AKZENTA INTERNATIONAL SA

Via Giuseppe Motta, 24
CH-6830 Chiasso

tel. +41 (0) 91 921 14 92
fax +41 (0) 91 682 85 20

e-mail: business@akzenta.com
www.akzenta.com

USt-IdNr. DE815221533
UID Nr. CHE-337.693.284 IVA
Part. IVA IT09701720964
Company no. CHE-337.693.284
Commercial Register: Switzerland, Canton Tessin

This catalogue was last updated in April 2018. For the most up to date version please go to akzenta.com/catalogue

<u>ABOUT US</u>	<u>4</u>	<u>DENTAL ROLLS</u>	<u>37</u>
<u>WHAT’S NEW 2018</u>	<u>6</u>	<u>HEADREST COVERS</u>	<u>38</u>
<u>PRODUCT RANGE EXPLAINED</u>	<u>8</u>	<u>SALIVA EJECTORS</u>	<u>41</u>
<u>APPLICATORS</u>	<u>10</u>	<u>MIXING PADS</u>	<u>42</u>
<u>APRONS</u>	<u>12</u>	<u>PATIENT TOWELS</u>	<u>44</u>
<u>ASPIRATOR TIPS</u>	<u>13</u>	<u>PROTECTIVE SLEEVES</u>	<u>50</u>
<u>BIB HOLDER</u>	<u>14</u>	<u>STERILISATION</u>	<u>52</u>
<u>CHAIR COVER</u>	<u>15</u>	<u>TISSUE NAPKINS</u>	<u>54</u>
<u>COOLING</u>	<u>16</u>	<u>TOOTHBRUSHES</u>	<u>55</u>
<u>CUPS</u>	<u>17</u>	<u>TRAY PAPER INSERTS</u>	<u>56</u>
<u>FACE MASKS</u>	<u>22</u>	<u>WET CREPE NAPKINS</u>	<u>57</u>
<u>GLOVES</u>	<u>26</u>	<u>PRODUCT SETS</u>	<u>58</u>
<u>AKZENTA STATISTICS</u>	<u>34</u>	<u>PRODUCT LINES IN DETAIL</u>	<u>59</u>
<u>HANDTOWELS</u>	<u>36</u>	<u>GET TO KNOW US</u>	<u>70</u>




ABOUT US

We at Akzenta care deeply about a perfect patient experience. Being on the market for almost 40 years, our partners trust us thanks to our values reliability, high quality, unique variety and fair prices. Akzenta has been manufacturing private label products for well-known international dental supply companies for over 30 years. We are a main contributor to a friendlier, more colourful dental practice atmosphere and are known trend-setters.



Akzenta is a trusted brand all over the world.



Private Label Service

We have been producing private label products for the EU's leading dental groups for almost 40 years. Market leaders like: Lifco, Plandent or Henry Schein.



Highly Trusted Brand

Customers across more than 30 countries trust in our cutting edge products and innovative dental articles. Our customer satisfaction is unrivaled in the industry.



Always One Step Ahead

We are trendsetters: whether it's our signature style of COCO Color Collection products, our beautiful new Paradise Line or our innovative Smily Series for kids: with Akzenta you're always on top of the market.



Highest Quality & Reliability

Akzenta provides you with a wide range of high quality products. Throughout the years Akzenta has established strong bonds with manufacturers and logistic partners to provide you with a smooth and consistent customer experience.

WHAT'S NEW

FOG STOP SWISS EDITION

No more foggy glasses! With the introduction of our brand new Fog Stop Flap you can concentrate on your work with crystal clear glasses.



BARRIER ENVELOPES

Our Barrier Envelopes protect your investment in digital imaging technology and are made to maximize patient comfort.

HYGIENIC SLEEVES

Convey supreme hygiene to your customers by using Hygienic Sleeves. Made from soft, but strong plastic, the hygienic sleeves work with a variety of holders.



FACEMASK BLACK

Try something new and exciting! Our Black Series was specifically made for professionals dealing with visible residue to convey cleanliness to their customers.



PARADISE EDITION

The Paradise Edition was made to create a cheerful and relaxed feeling in your dental practice. For happier patients – available in 4 colors.



STYLE NITRILE GLOVE LAGUNA BLUE

Our popular STYLE NITRILE Series just got supplemented with a new and great color: Laguna Blue. This glove will fit perfectly into the discreet look of a modern practice.



THE AKZENTA PRODUCT RANGE EXPLAINED

Top Selection

The Top Selection is Akzenta's high quality line aimed at giving you superb disposables at very competitive prices. Experience quality and comfort in a new dimension.



Swiss Edition

Swiss Quality is recognised for its precision and dedication to excellence worldwide. Delight your customer with our premium Swiss Edition products.



Top Selection+

The Top Selection+ delivers the most sophisticated selection of disposables. Aimed at those who demand nothing but the best.



Paradise Edition

Create a more pleasant and relaxed atmosphere with our Paradise Edition. The friendly motives bring a sense of vacation to your dental practice.



Smily Series

The Smily Series aims at making visits to the dentist's office more joyful for children. Put your younger patients at ease with our cute designs.



Black Series

Be individual and make a visual statement with our Black Series. The pitch black color is both practical and appealing.



COCO Color Collection

The COCO Color Collection makes dental practices more colourful and your patients smile. Our 11 most popular products in 7 vibrant colors.



Pure Edition

Our bestselling Pure Edition Gloves are the softest gloves you will find. Its silky smooth surface makes long wearing periods a breeze.



APPLICATORS

3 NEW
COLORS



Top Selection+



TOP STICKS PLUS MICRO APPLICATORS | REGULAR



Med Blue
MA-AZ-MB



Burgundy
MA-AZ-WR



Fresh Green
MA-AZ-FG



Fuchsia
MA-AZ-F



Lilac
MA-AZ-L



Orange
MA-AZ-O



Pink
MA-AZ-P

| FINE



Yellow
MA-AZ-Y

| SUPER FINE



White
MA-AZ-W

Strong and non-absorbent, non-linting applicators, with bendable neck and flexible plastic handle for pinpoint accuracy application of any liquids.

PCS Dispenser	Dispenser Box	Boxes Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Pallet Size CM	Pallet Weight KG
100	4	40	52 x 42 x 24.5	15.5	33	80 x 120 x 245	531.5



APRONS

TOP SAFE PLUS APRONS | REGULAR



Fuchsia
PATU-AZ-F



White
PATU-AZ-W



Yellow
PATU-AZ-Y



Peach
PATU-AZ-PE



Orange
PATU-AZ-O



Pink
PATU-AZ-P



Lilac
PATU-AZ-L



Burgundy
PATU-AZ-WR



Med Blue
PATU-AZ-MB



Fresh Green
PATU-AZ-FG



Light Blue
PATU-AZ-B



Green
PATU-AZ-M

Strong protection against contamination, highly absorbent on the surface. Protective polythene backing 1-ply tissue and 1-ply PE-foil.


Size: Regular | 50 x 60 cm

PCS Roll	Rolls Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layer Pallet	Pallet Size CM	Pallet Weight KG
80	6	35 x 23.5 x 26.2	7.8	60	10	6	80 x 120 x 180	483

TOP SAFE PLUS | APRONS WITH POCKET | REGULAR



Blue 38x68
BIB3L-AZ-B




White 38x68
BIB3L-AZ-W



Blue 38x48
BIB3M-AZ-B



White 38x48
BIB3M-AZ-W



ASPIRATOR TIPS

TOP TIPS & TOP TIPS | KIDS



White
ABSA-AZ-W
Kids:
ABSB1-AZ-W



Yellow
ABSA-AZ-Y
Kids:
ABSB1-AZ-Y



Orange
ABSA-AZ-O
Kids:
ABSB1-AZ-O



Red
ABSA-AZ-R
Kids:
ABSB1-AZ-R



Burgundy
ABSA-AZ-WR
Kids:
ABSB1-AZ-WR



Pink
ABSA-AZ-P
Kids:
ABSB1-AZ-P



Lilac
ABSA-AZ-L
Kids:
ABSB1-AZ-L

Synthetic material, high production standard, non-slip surface, low noise, well formed, therefore easy and perfect handling. Matt-finished surface (no interference from light reflections) autoclavable up to 134°C CE

Consumable units:
10 Pieces | Polybag

TOP TIPS Size:
Length 124 mm | Diameter 16 mm | Polybag
Size 170x125x017 mm | Polybag Weight 0.08 Kg

TOP TIPS KIDS Size:
Length 108 mm | Diameter 16 mm | Polybag
Weight 0.06 Kg | Polybag Size 145x110x015 mm



Fuchsia
ABSA-AZ-F
Kids:
ABSB1-AZ-F



Med Blue
ABSA-AZ-MB
Kids:
ABSB1-AZ-MB



Blue
ABSA-AZ-B
Kids:
ABSB1-AZ-B



Fresh Green
ABSA-AZ-FG
Kids:
ABSB1-AZ-FG



Green
ABSA-AZ-G
Kids:
ABSB1-AZ-G



BIB HOLDER

TOP EASY





BH-N

TOP EASY PLUS



BH-AZ

The new hygienic, flexible and disposable bib holder solution.



PCS Dispenser	Dispenser Carton	Boxes Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	PCS Carton	Layer Pallet	Pallet Size CM	Pallet Weight KG
250	25	30	48 x 32 x 33	5	35	5	6 250	7	80 x 120 x 240	135

CHAIR COVER

TOP COVER PLUS

- Protects your dental chair from dirt and deterioration
- Comes in 2 or 3 pcs set
- Fits in with our coco color collection

AVAILABLE COLORS



BLUE



GREEN



FUCHSIA



FRESH GREEN



LILAC



MED BLUE



ORANGE



YELLOW



3PCS

1. Headrest Cover
2. Back
3. Leg

CC3-AZ-B



2PCS

2. Back
3. Leg

CC2-AZ-B



Article number	Color	PCS Unit	Units Carton	Carton Size	Carton Weight Kg	Cartons Layer	Layers Pallet	Cartons Pallet	Pallet Size CM	Pallet Weight Kg
CC3-AZ-B	Blue	25 x 3	4 x (25 x 3)	58 x 40 x 32	8	4	6	24	80 x 120 x 207	212
CC3-AZ-G	Green	25 x 3	4 x (25 x 3)	58 x 40 x 32	8	4	6	24	80 x 120 x 207	212
CC3-AZ-F	Fuchsia	25 x 3	4 x (25 x 3)	58 x 40 x 32	8	4	6	24	80 x 120 x 207	212
CC3-AZ-FG	Fresh Green	25 x 3	4 x (25 x 3)	58 x 40 x 32	8	4	6	24	80 x 120 x 207	212
CC3-AZ-L	Lilac	25 x 3	4 x (25 x 3)	58 x 40 x 32	8	4	6	24	80 x 120 x 207	212
CC3-AZ-MB	Med Blue	25 x 3	4 x (25 x 3)	58 x 40 x 32	8	4	6	24	80 x 120 x 207	212
CC3-AZ-O	Orange	25 x 3	4 x (25 x 3)	58 x 40 x 32	8	4	6	24	80 x 120 x 207	212
CC3-AZ-Y	Yellow	25 x 3	4 x (25 x 3)	58 x 40 x 32	8	4	6	24	80 x 120 x 207	212
CC2-AZ-B	Blue	25 x 2	4 x (25 x 2)	58 x 40 x 32	5.5	5	6	30	80 x 120 x 195	185
CC2-AZ-G	Green	25 x 2	4 x (25 x 2)	58 x 40 x 32	5.5	5	6	30	80 x 120 x 195	185

COOLING

Content ml	Cans Carton	Carton Size CM	Carton Weight Kg	Layers Carton	Layers Pallet	Carton Pallet	Pallet Size CM	Pallet Weight KG
200	12	21 x 16 x 19.5	2	21	8	168	80x120x171	440



400-ICET-L02
100-ICET-L02

By pressing and shaking the **TOP FROST bag** an internal endothermic reaction between the water and the azotic salt compound is triggered. As a result, the temperature drops down to -5°C.

External bag made of non-woven fabric with a triple heat sealing to guarantee highest resistance, cold storage and impermeability. The non-woven fabric makes sure the temperature is evenly distributed. Patent protected.

Mixture of azotic salt according to EU-norm CE 1907/2006.


Article Number	Format cm	Bag Carton	Carton Size CM	Carton Weight Kg	Cartons Layers	Layers Pallet	Carton Pallet	Pallet Size CM	Pallet Weight KG
400-ICET-L02	13.5 x 18	10	29 x 18.5 x 15	3	16	8	128	80 x 120 x 135	404
100-ICET-L02	13.5 x 18	25	38 x 25 x 15.6	7.3	9	9	81	80 x 120 x 155	611.3



COOL-AZ

TOP COOL Spray cools up to -45 °C and vaporizes immediately, residue-free. It is ideal to instantly freeze foam pellets and dental rolls to perform an indirect vitality check on the tooth. It is not meant for direct oral use.

It can also be used for rapid cooling of impressions and thermoforming films (medical device class 1). A precise application tube serves for pinpoint application. Top Cool Spray comes in peppermint flavor, filled with propellant in cosmetic quality.



CUPS



Our high quality dental cups are made of food safe polypropylene or polystyrene. No crunching and comfortably smooth – but very resistant and solid – are qualities we guarantee. The cups are colorful and shiny.

We use the latest technology and most modern production without any human contact to ensure absolute hygiene.

Suitable for all standardized cup dispensers.
100% recyclable
Burning without any residue.

CUPS

TOP CUPS WHITE PP



Article Number	Content	PCS Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Carton Layers	Layer Pallet	Pallet Size CM	Pallet Weight KG
CUPP150-1-N-W	150 ml	1000	39.5 x 28 x 31.5	3	63	9	7	80 x 120 x 236	209
CUPP150-3-N-W	150 ml	3000	38 x 38 x 51.5	8.5	24	6	4	80 x 120 x 221	228
CUPP150-3-N-W	210 ml	1000	39.5 x 28 x 31.5	3	63	9	7	80 x 120 x 236	209
CUPP200LL-15-N-W	210 ml	1500	38 x 38 x 34	4.5	36	6	6	80 x 120 x 220	195
CUPP200LL-3-N-W	210 ml	3000	38 x 38 x 30	8.5	24	6	4	80 x 120 x 221	228

TOP CUPS WHITE PS



Article Number	Content	PCS Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Carton Layers	Layer Pallet	Pallet Size CM	Pallet Weight KG
108-CUPS-656	180 ml	2000	39 x 29 x 59	7	36	8	4	80 x 120 x 250	272
CUPS200-2-N-W	200 ml	2000	39 x 29 x 59	7.3	32	8	4	80 x 120 x 193	254
CUPS225-2-N-W	225 ml	2000	39 x 29 x 59	7.5	32	8	4	80 x 120 x 193	260



TOP CUPS PP



White
CUPP3-N-W
CUPP15-N-W
CUPP1-N-W



Yellow
CUPP3-N-Y
CUPP15-N-Y
CUPP1-N-Y



Orange
CUPP3-N-O
CUPP15-N-O
CUPP1-N-O



Lilac
CUPP3-N-L
CUPP15-N-L
CUPP1-N-L



Viola
CUPP3-N-VI
CUPP15-N-VI
CUPP1-N-VI



Burgundy
CUPP3-N-WR
CUPP15-N-WR
CUPP1-N-WR



Fresh Green
CUPP3-N-FG
CUPP15-N-FG
CUPP1-N-FG



Pink
CUPP3-N-P
CUPP15-N-P
CUPP1-N-P



Light Blue
CUPP3-N-LB
CUPP15-N-LB
CUPP1-N-LB



Dark Green
CUPP3-N-DG
CUPP15-N-DG
CUPP1-N-DG



Fuchsia
CUPP3-N-F
CUPP15-N-F
CUPP1-N-F



Dark Blue
CUPP3-N-DB
CUPP15-N-DB
CUPP1-N-DB



Transparent
CUPP3-N-TR
CUPP15-N-TR
CUPP1-N-TR



Med Blue
CUPP3-N-MB
CUPP15-N-MB
CUPP1-N-MB



Red
CUPP3-N-R
CUPP15-N-R
CUPP1-N-R



Mint Green
CUPP3-N-MG
CUPP15-N-MG
CUPP1-N-MG

PCS Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Carton Layers	Layer Pallet	Pallet Size CM	Pallet Weight KG
3000	38 x 38 x 51.5	8.5	24	6	4	80 x 120 x 221	228
1500	38 x 38 x 34	4.5	36	6	6	80 x 120 x 220	195
1000	38 x 26.5 x 30	3	63	9	7	80 x 120 x 225	209



PCS Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Carton Layers	Layer Pallet	Pallet Size CM	Pallet Weight KG
1000	38 x 26.5 x 30	3	63	9	7	80 x 120 x 225	209

COCO Color Collection

STYLE CUPS BICOLOR



PCS Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
1200	40 x 40 x 32	4.6	36	6	6	80 x 120 x 207	187



SMILY CUPS



White
CUPP1-SMILE

The high quality Smily Cup made of food safe, authentic polypropylene assures: no crunching, comfortably smooth, but very resistant and solid when filled with water.

PCS Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
1000	38 x 26.5 x 30	3	63	9	7	80 x 120 x 225	209



FACE MASKS



FACE MASKS

TOP MASK TYPE IIR 3-PLY FACE MASK | NON LATEX



White
FMEB2AZ-W



Yellow
FMEB2AZ-Y



Orange
FMEB2AZ-O



Pink
FMEB2AZ-P



Lilac
FMEB2AZ-L



Green
FMEB2AZ-G



Fresh Green
FMEB2AZ-FG



Fuchsia
FMEB2AZ-F



Blue
FMEB2AZ-B



Burgundy
FMEB2AZ-WR



Med Blue
FMEB2AZ-MB

Our Top Mask is distinguished by its optimum fit. It provides a high wearing comfort even during extended wearing time. Only eudermic materials are used. Filtration efficiency $\geq 98\%$ acc. EN14683 type II R

Type	PCS Dispenser	DIspensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
TOP MASK	50	20	52 x 19 x 34	4.5	40	8	5	80 x 120 x 185	185

TOP MASK PLUS | COOLBREATH - TYPE IIR



White
FMEBC2+AZ-W



Yellow
FMEBC2+AZ-Y



Orange
FMEBC2+AZ-O



Pink
FMEBC2+AZ-P

This face mask offers enhanced performance against micro-dust, anti-spray and viral particles. It provides a high wearing comfort even during extended wearing time. Only eudermic materials are used. Filtration efficiency $\geq 98\%$ acc. EN14683. IIR – resistant to drops, very strong and high quality.



Lilac
FMEBC2+AZ-L



Green
FMEBC2+AZ-G



Fresh Green
FMEBC2+AZ-FG



Blue
FMEBC2+AZ-B

PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
50	20	52 x 20 x 38	4.6	40	8	5	80 x 120 x 205	204



Black
FMEB2AZ-BK

The Face Mask Black Series is our newest development for progressive and state-of-the-art practices. Only eudermic materials are used. Filtration efficiency $\geq 98\%$ acc. EN14683 type II R



White
FMEB2+F-AZ-W



Green
FMEB2+F-AZ-G



Blue
FMEB2+F-AZ-B

Our Facemask Sensitive Fog Stop Swiss Edition is distinguished by its revolutionary flap, which prevents air from making your glasses foggy while working. This makes it one of the most modern and capable innovations we recently designed.

Smily Series

Paradise Edition



White
FMEB2+SMILY

The SMILY Mask is designed specifically for younger patients and goes well with our other SMILY products. Only eudermic materials are used. Filtration efficiency $\geq 98\%$ acc. EN14683 type II R



Fuchsia
FMEB2AZ-PAR-F



Fresh Green
FMEB2AZ-PAR-FG



Lilac
FMEB2AZ-PAR-L



Light Blue
FMEB2AZ-PAR-B

Bring a sense of vacation to your dental practice. The new Paradise Edition sets a new standard in patient comfort and makes visits to the dentist a breeze.

PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
50	20	52 x 19 x 34	4.5	40	8	5	80 x 120 x 185	200



GLOVES



Top Selection

TOP TOUCH | LATEX | POWDERED



Natural White

Article Number	Size	PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
LATX-1-AZ-XS	XS	100	10	38 x 25 x 25	5.6	72	9	7	80 x 120 x 190	372
LATX-1-AZ-S	S	100	10	38 x 25 x 25	6.0	72	9	7	80 x 120 x 190	398
LATX-1-AZ-M	M	100	10	38 x 25 x 25	6.4	72	9	7	80 x 120 x 190	423
LATX-1-AZ-L	L	100	10	38 x 25 x 25	6.8	72	9	7	80 x 120 x 190	448

This powdered glove has a very good fit, is soft and easy to put on and to take off. Its textured surface ensures an excellent grip and a better sense of touch. TOP TOUCH Latex Powdered is made of hygienic raw material, free of thiuram, plasticier, animal products and any other biological substances. It is biodegradable and environmentally friendly. It is supplied in an attractive packaging.

TOP TOUCH | LATEX | EXTRA SENSITIVE



White

Article Number	Size	PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
LATX7-AZ-XS	XS	100	10	25 x 34 x 25	5.6	63	9	7	80 x 120 x 190	372
LATX7-AZ-S	S	100	10	25 x 34 x 25	6	63	9	7	80 x 120 x 190	398
LATX7-AZ-M	M	100	10	25 x 34 x 25	6.4	63	9	7	80 x 120 x 190	423
LATX7-AZ-L	L	100	10	25 x 34 x 25	6.8	63	9	7	80 x 120 x 190	448

The TOP TOUCH Latex Extra Sensitive is a glove developed to be easy to put on, has a good fit and is extra flexible. Its surface is smooth. The TOP TOUCH is powder-free and polymer-coated. The protein content is less than 50 µg/g. It is made of hygienic raw material, free of thiuram, plasticier, animal products and any other biological substances. It is biodegradable and environmentally friendly.

GLOVES



White

TOP TOUCH PLUS | NITRILE ULTRASOFT

Article Number	Size	PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
NITR2S-AZ-W-XS	XS	100	10	27 x 25.5 x 24.5	3.5	84	12	7	80 x 120 x 187	314
NITR2S-AZ-W-S	S	100	10	27 x 25.5 x 24.5	3.7	84	12	7	80 x 120 x 187	330.8
NITR2S-AZ-W-M	M	100	10	27 x 25.5 x 24.5	3.9	84	12	7	80 x 120 x 187	347.6
NITR2S-AZ-W-L	L	100	10	27 x 25.5 x 24.5	4.1	84	12	7	80 x 120 x 187	364.4



Lavender

TOP TOUCH PLUS | NITRILE LAVENDER

Article Number	Size	PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
NITR2S-AZ-LA-XS	XS	100	10	31 x 26 x 24.5	4	84	12	7	80 x 120 x 187	356
NITR2S-AZ-LA-S	S	100	10	31 x 26 x 24.5	4.3	84	12	7	80 x 120 x 187	381
NITR2S-AZ-LA-M	M	100	10	31 x 26 x 24.5	4.5	84	12	7	80 x 120 x 187	398
NITR2S-AZ-LA-L	L	100	10	31 x 26 x 24.5	4.8	84	12	7	80 x 120 x 187	423



Magic Blue

TOP TOUCH PLUS | NITRILE | MAGIC BLUE

Article Number	Size	PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
NITR2S-AZ-MBL-XS	XS	100	10	31 x 26 x 24.5	4	84	12	7	80 x 120 x 187	356
NITR2S-AZ-MBL-S	S	100	10	31 x 26 x 24.5	4.3	84	12	7	80 x 120 x 187	381
NITR2S-AZ-MBL-M	M	100	10	31 x 26 x 24.5	4.5	84	12	7	80 x 120 x 187	398
NITR2S-AZ-MBL-L	L	100	10	31 x 26 x 24.5	4.8	84	12	7	80 x 120 x 187	423

Due to its thinness, a high wearing comfort is guaranteed. The glove is also very flexible and has a natural feeling. Best grip is guaranteed through its roughened fingertips.



Light Beige

TOP TOUCH | LATEX | SWISS EDITION

Article Number	Size	PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
LATX3-XS	XS	100	10	34 x 25 x 25	6	63	9	7	80 x 120 x 193	398
LATX3-S	S	100	10	34 x 25 x 25	6.5	63	9	7	80 x 120 x 193	429
LATX3-M	M	100	10	34 x 25 x 25	7	63	9	7	80 x 120 x 193	461
LATX3-L	L	100	10	34 x 25 x 25	7.5	63	9	7	80 x 120 x 193	492
LATX3-XL	XL	100	10	34 x 25 x 25	8	63	9	7	80 x 120 x 195	524

This powder-free glove has a very good fit, is soft and easy to put on and to take off. Its textured surface ensures an excellent grip and a better sense of touch. The TOP TOUCH LATEX SWISS EDITION is made of hygienic raw material, free of thiuram, plasticier, animal tissue and any other biological substances. It is biodegradable and environmentally friendly. It is supplied in an attractive packaging.



GLOVES

STYLE SERIES | NITRILE



Laguna Blue
STYLE1-AZ-LB-XS
STYLE1-AZ-LB-S
STYLE1-AZ-LB-M
STYLE1-AZ-LB-L

The STYLE NITRILE is a super sensitive glove in fashionable colors. It provides a great comfort and excellent fit. The STYLE NITRILE is extra thin and ultrasensitive, soft, smooth and easy to put on and take off. Its micro-textured fingertips ensure an excellent tactile sensitivity.

Size	PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
XS	100	10	26 x 25.5 x 25.2	3.5	84	12	7	80 x 120 x 194	314
S	100	10	26 x 25.5 x 25.2	3.8	84	12	7	80 x 120 x 194	339
M	100	10	26 x 25.5 x 25.2	4	84	12	7	80 x 120 x 194	356
L	100	10	26 x 25.5 x 25.2	4.5	84	12	7	80 x 120 x 194	398



Pink
STYLE1-AZ-P-XS
STYLE1-AZ-P-S
STYLE1-AZ-P-M
STYLE1-AZ-P-L



Med Blue
STYLE1-AZ-MB-XS
STYLE1-AZ-MB-S
STYLE1-AZ-MB-M
STYLE1-AZ-MB-L



Fresh Green
STYLE1-AZ-FG-XS
STYLE1-AZ-FG-S
STYLE1-AZ-FG-M
STYLE1-AZ-FG-L



Fuchsia
STYLE1-AZ-F-XS
STYLE1-AZ-F-S
STYLE1-AZ-F-M
STYLE1-AZ-F-L



Orange
STYLE1-AZ-O-XS
STYLE1-AZ-O-S
STYLE1-AZ-O-M
STYLE1-AZ-O-L

STYLE SERIES | LATEX



Pink
STYLE-AZ-P-XS
STYLE-AZ-P-S
STYLE-AZ-P-M
STYLE-AZ-P-L



Med Blue
STYLE-AZ-B-XS
STYLE-AZ-B-S
STYLE-AZ-B-M
STYLE-AZ-B-L



Fresh Green
STYLE-AZ-G-XS
STYLE-AZ-G-S
STYLE-AZ-G-M
STYLE-AZ-G-L



Fuchsia
STYLE-AZ-F-XS
STYLE-AZ-F-S
STYLE-AZ-F-M
STYLE-AZ-F-L



Orange
STYLE-AZ-O-XS
STYLE-AZ-O-S
STYLE-AZ-O-M
STYLE-AZ-O-L



Lilac
STYLE-AZ-L-XS
STYLE-AZ-L-S
STYLE-AZ-L-M
STYLE-AZ-L-L

Size	PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
XS	100	10	39.5 x 26 x 25.5	6	63	9	7	80 x 120 x 194	398
S	100	10	39.5 x 26 x 25.5	6.5	63	9	7	80 x 120 x 194	430
M	100	10	39.5 x 26 x 25.5	7.5	63	9	7	80 x 120 x 194	492

Our STYLE Gloves made of premium super elastic latex are powderfree and ambidextrous. Available in 6 attractive colors. Provides excellent grip.



GLOVES



Black

TOP TOUCH PLUS | NITRILE BLACK

Article Number	Size	PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
NITR2-AZ-BK-S	S	100	10	26 x 26 x 24.5	3.6	84	12	7	80 x 120 x 186	322
NITR2-AZ-BK-M	M	100	10	26 x 26 x 24.5	4.73	84	12	7	80 x 120 x 186	417
NITR2-AZ-BK-L	L	100	10	26 x 26 x 24.5	5.23	84	12	7	80 x 120 x 186	459
NITR2-AZ-BK-XL	XL	100	10	26 x 26 x 24.5	5.53	84	12	7	80 x 120 x 186	484

Due to its thinness, a high wearing comfort is guaranteed. The glove is also very flexible and has a natural feeling. Best grip is guaranteed by its roughened fingertips.



Black

TOP TOUCH PLUS | LATEX BLACK

Article Number	Size	PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
LATX-2-AZ-BK-S	S	100	10	39.5 x 26 x 25.5	6.5	63	9	7	80 x 120 x 193	429
LATX-2-AZ-BK-M	M	100	10	39.5 x 26 x 25.5	7	63	9	7	80 x 120 x 193	461
LATX-2-AZ-BK-L	L	100	10	39.5 x 26 x 25.5	7.5	63	9	7	80 x 120 x 193	492
LATX-2-AZ-BK-XL	XL	100	10	39.5 x 26 x 25.5	8	63	9	7	80 x 120 x 193	524

The TOP TOUCH PLUS Latex Black is easy to put on, has a good fit and high wearing comfort is guaranteed due to its fineness. TOP TOUCH PLUS Latex Black is powder-free and skin friendly. Its fingertips are textured and enable excellent instrument handling. The TOP TOUCH Latex Black is very elastic and has a natural touch.



White

PURE EDITION | NITRILE SOFT

Article Number	Size	PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
NITR-3-09XS	XS	100	10	33.5 x 25.2 x 25	4.2	63	9	7	80 x 120 x 190	284
NITR-3-09S	S	100	10	33.5 x 25.2 x 25	4.6	63	9	7	80 x 120 x 190	310
NITR-3-09M	M	100	10	33.5 x 25.2 x 25	5	63	9	7	80 x 120 x 190	335
NITR-3-09L	L	100	10	33.5 x 25.2 x 25	5.4	63	9	7	80 x 120 x 190	360

This glove is ambidextrous and has a very good fit. The fingers are textured. The glove provides great tactile sensitivity and an excellent grip. Non-sterile. Attractive packaging.




White

PURE EDITION | LATEX SENSITIVE

Article Number	Size	PCS Dispenser	Dispensers Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
LATX-3-19XS	XS	100	10	34 x 25 x 25	6	63	9	7	80 x 120 x 190	398
LATX-3-19S	S	100	10	34 x 25 x 25	6.5	63	9	7	80 x 120 x 190	429
LATX-3-19M	M	100	10	34 x 25 x 25	7	63	9	9	80 x 120 x 190	461
LATX-3-19L	L	100	10	34 x 25 x 25	7.5	63	9	9	80 x 120 x 190	492

This powder-free glove has a very good fit, is soft and easy to put on and to take off. Its textured surface ensures an excellent grip and a better sense of touch. The PURE LATEX SENSITIVE is made of hygienic raw material, free of thiuram, plasticier, animal tissue and any other biological substances. It is biodegradable and environmentally friendly after use. It is supplied in an attractive packaging.

AKZENTA STATISTICS




Founded in 1979




**Market leader for
Tray Paper Inserts
and Instrument
Table Linings**

**Trusted quality from Switzerland
Proudly based in Tessin, Switzerland since 1989**




Sterilisation Paper made in Switzerland



**Tray Filter paper:
100% biodegradable**

12 Mio. € revenue each year, 14 employees





**Over 1.000 tons
of Tray Paper sold
every year**



**More than 200 customers
in over 30 countries**

**Over 200 Mio. cups
sold to customers
every year**



**Only company worldwide to produce
a whole product line in the same exact
pantone color**

**Big selection of products on stock
that are shipping within 24 hours in Germany
/within 4 days in Europe**



HANDTOWELS


TOP TOWELS PLUS, HANDFOLD TOWELS Z-FOLD



White

Folded hand towels
Recycled Paper Mix
Double layered
Extra soft tissue
Highly absorbent for efficient hand drying
Suitable for all Standard Dispensers
HAT- Highly Absorbing Tissue

Article Number	Size	PCS Pack	Pack Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
FAHA-AZ-ZW2Z	24 x 21	160	20	60 x 40 x 25	6.8	36	4	9	80 x 120 x 240	264




DENTAL ROLLS

TOP ROLLS PLUS



Our Top Rolls plus are double protected for best hygienic standards in a special open and close zip bag. Each dispenser consists of 20 bundles of 300gr. The Rolls have a width off 12, 10 or 8mm and a length of 40mm.

Article Number	Size	Width mm	Length mm	Bundle Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
WARO+-AZ-3	3	12	40	20 x 300g	40 x 40 x 35	6.5	24	6	4	80 x 120 x 155	176
WARO+-AZ-2	2	10	40	20 x 300g	40 x 40 x 35	6.5	24	6	4	80 x 120 x 155	176
WARO+-AZ-1	1	8	40	20 x 300g	40 x 40 x 35	6.5	24	6	4	80 x 120 x 155	176



HEADREST COVERS



Top Selection



TOP COVER TISSUE-PE HEADREST COVER



White
KOPF-0-33W



Yellow
KOPF-0-33Y



Rosa
KOPF-0-33P



Peach
KOPF-0-33PE



Lilac
KOPF-0-33L



Mint Green
KOPF-0-33ME

Absorbing (on top) and protecting (back side). 1-ply extra-strong (26g/m²), highly absorbent tissue of pure cellulose and 1-ply PE-foil, latexfree. Tissue 26 g/m2 + polyethylene.



Fresh Green
KOPF-0-33FG



Light Blue
KOPF-0-33LB



Med Blue
KOPF-0-33MB



Orange
KOPF-0-33O



Fuchsia
KOPF-0-33F



Burgundy
KOPF-0-33WR

Pieces Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
500	36 x 27 x 19	3.9	90	9	10	80 x 120 x 205	371

HEADREST COVERS

ECO SILK PAPER HEADREST COVER



White
KOPF28C-AZ
KOPF36C-AZ

Made out of Tissue Paper 25g/m2
our best price/performance-ratio!

Size	Pieces Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Layers Pallet	First 8 Layers: Cartons Layer	Last 9 th Layer: Cartons Layer	Pallet Size CM	Pallet Weight KG
28 x 28	1,000	29 x 29 x 9	3.6	100	9	12	4	80 x 120 x 96	380
28 x 36	1,000	37 x 29 x 9	4.2	100	9	12	4	80 x 120 x 96	440



TOP COVER STERILISATION-PAPER HEADREST COVER

The hygienic solution
Environmentally friendly production without use of any harmful
materials



White
HCSP28-AZ-W
HCSP36-AZ-W



Yellow
HCSP28-AZ-Y
HCSP36-AZ-Y



Green
HCSP28-AZ-G
HCSP36-AZ-G



Blue
HCSP28-AZ-B
HCSP36-AZ-B

Size	Pieces Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
28 x 28	200	37.6 x 29.6 x 10.2	1.16	120	8	15	80 x 120 x 168	155
28 x 36	200	37.6 x 29.6 x 10.2	1.44	120	8	15	80 x 120 x 168	188

SALIVA EJECTORS



White
SE2-AZ-W
SE-AZ-W



Transparent
SE2-AZ-TR
SE-AZ-TR



Yellow
SE2-AZ-Y
SE-AZ-Y



Orange
SE2-AZ-O
SE-AZ-O



Light Pink
SE2-AZ-P
SE-AZ-P



Fuchsia
SE2-AZ-F
SE-AZ-F



Red
SE2-AZ-R
SE-AZ-R

Medical PVC with fixed weld copper- or zinc-wire
based on splendid adaptation of the single components. Easily
bendable and firm, clean cut edge, strong draw tip with sidewise
deairing for avoiding suction on the mucus membrane. CE.



Lilac
SE2-AZ-L
SE-AZ-L



Burgundy
SE2-AZ-WR
SE-AZ-WR



Fresh Green
SE2-AZ-FG
SE-AZ-FG



Med Blue
SE2-AZ-MB
SE-AZ-MB



Dark Green
SE2-AZ-DG
SE-AZ-DG



Dark Blue
SE2-AZ-DB
SE-AZ-DB



Light Blue
SE2-AZ-LB
SE2-AZ-LB



Mint Green
SE2-AZ-ME
SE-AZ-ME

Format mm	PCS Dispenser	PCS Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
125 (SE2)	100	10	29 x 18 x 20.5	4	144	16	9	80 x 120 x 200	596
145 (SE)	100	10	29 x 18 x 20.5	4.6	144	16	9	80 x 120 x 200	682



MIXING PADS



TOP MIX MIXING PADS PARCHMENT ANTI-SLIDE

Splendidly suitable for mixing of cements and casting materials
Very reasonably priced in 3-block version
Fixed ground during mixing process by anti-slide / non slip foam base

Article Number	Format mm	Sheets Block	Blocks Carton	Carton Size CM	Carton Weight Kg
MAPG-6262	62 x 62	50	275	31 x 21.5 x 11	4
MAPG-70100	70 x 100	75	120	31 x 21.5 x 11	4.35

TOP MIX MIXING PADS PE-COATED PAPER ANTI-SLIDE

Splendidly suitable for mixing of cements and casting materials
Very manageable reasonably priced in 3-block version
Fixed ground during mixing process by anti-slide / non slip foam base

Article Number	Format mm	Papers Block	Blocks Carton	Carton Size CM	Carton Weight Kg
MIPA-7095	70 x 95	50	81	31 x 22.5 x 11	6.79
MIPA-125190	125 x 190	25	48	31 x 22.5 x 11	6.79

TOP MIX MIXING PADS PARCHMENT


Splendidly suitable for mixing of cements and casting materials
Very manageable and reasonably priced in 4-block version

Article Number	Format mm	Sheets Block	Blocks Carton	Carton Size CM	Carton Weight Kg
MIPG-70100	70 x 100	75	120	31 x 21.5 x 11	4.35
MIPG-K4050	40 x 40	75 x 3	167	31 x 22.5 x 11	4.65
MIPG-K70100	70 x 100	75 x 3	40	31 x 21.5 x 11	4.35

TOP MIX MIXING PADS PE-COATED PAPER

Splendidly suitable for mixing of cements and casting materials.
Very manageable reasonably priced in 3-block version

Article Number	Format mm	Sheets Block	Blocks Carton	Carton Size CM	Carton Weight Kg
MIPE-K7095	70 x 95	50 x 3	27	31 x 22.5 x 11	6.79
MIPE-K125190	125 x 190	50 x 3	24	31 x 22.5 x 11	7.23
MIPE-K150240	150 x 240	50 x 3	12	31 x 22.5 x 11	5.6



PATIENT TOWELS



Top Selection



TOP SAFE PATIENT TOWELS 3-LAYERS



White
TOPA3+AZ-W



Yellow
TOPA3+AZ-Y



Peach
TOPA3+AZ-PE



Orange
TOPA3+AZ-O



Pink
TOPA3+AZ-P



Lilac
TOPA3+AZ-L



Fresh Green
TOPA3+AZ-FG



Light Blue
TOPA3+AZ-LB



Mint Green
TOPA3+AZ-G



Dark Blue
TOPA3+AZ-DB

Absorbing (on top) and protecting (back side). 2-ply extra-strong (17g/m2), highly absorbent tissue of pure cellulose and 1-ply PE-foil, latexfree.
2 x 17 g tissue + PE
The special features of the materials offer strong absorbtion and waterproofing to assure premium comfort of the patient.
Size: 33 x 45 cm



Pieces Pack	Pieces Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
50	10	33 x 23.8 x 28.5	3.5	70	10	7	80 x 120 x 215	265

PATIENT TOWELS

TOP SAFE PATIENT TOWELS 2-LAYERS



White
TOPA2+N-W

Yellow
TOPA2+N-Y

Peach
TOPA2+N-PE

Orange
TOPA2+N-O

Pink
TOPA2+N-P

Lilac
TOPA2+N-L

Burgundy
TOPA2+N-WR

The special features of the materials offer great absorption and waterproofing to assure premium comfort to the patient. Latex free, Size: 33 x 45 cm, 26 g/qm Tissue + 8 g/qm PE.



Fuchsia
TOPA2+N-F

Fresh Green
TOPA2+N-FG

Light Blue
TOPA2+N-B

Mint Green
TOPA2+N-ME

Med Blue
TOPA2+N-MB

Black
TOPA2+N-BK

Pieces Pack	Pieces Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
50	10 x 50	33 x 23.8 x 28.5	2.65	80	10	8	80 x 120 x 243	232

TOP SAFE PLUS



White
BAST37-AZ-W

Yellow
BAST37-AZ-Y

Orange
BAST37-AZ-O

Lilac
BAST37-AZ-L

Mint Green
BAST37-AZ-G

Med Blue
BAST37-AZ-B

Dark Blue
BAST37-AZ-DB




Fresh Green
BAST37-AZ-FG

Peach
BAST37-AZ-PE

Extra Soft, suitable as mouth wipe after its use. Layered with a PE Cover (Safe Back). Active protection from unwanted drops. Size: 37 x 37 cm, 26 g/qm Tissue + 8 g/qm PE

Pieces Pack	Pieces Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
50	500	37.5 x 37.5 x 13	3.5	96	6	16	80 x 120 x 223	356



SMILY TOWELS

Absorbing (on top) and protecting (back side).
1-ply extra-strong (26g/m2), highly absorbent tissue of pure cellulose and 1-ply PE-folie (8g/m2).
The special features of the materials offer high absorbtion and waterproofing to assure premium comfort to the patient.
Size: 33 x 37 cm.



White
BAST33-SMILY

Pieces Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
500	38.5 x 34.5 x 13	2.65	96	6	16	80 x 120 x 223	275

PARADISE TOWELS

Absorbing (on top) and protecting (back side). 1-ply extra-strong (26g/m2), highly absorbent tissue of pure cellulose and 1-ply PE-foil (8g/m2).
The special features of the materials offer high absorption and waterproofing to assure premium comfort to the patient.
Size: 33 x 37 cm.



Fuchsia
BAST33-PAR-F



Fresh green
BAST33-PAR-FG



Lilac
BAST33-PAR-L



Light Blue
BAST33-PAR-B

Pieces Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight Kg
500	38.5 x 34.5 x 13	2.65	96	6	16	80 x 120 x 223	275



PROTECTIVE SLEEVES



Top Selection+



BARRIER ENVELOPES



Size 0



Size 1



Size 2



Size 3

Protect your investment in digital imaging technology with our unique barrier envelopes, designed to fit all Dentsply, Kodak, Sirona and Air Techniques phosphor plates. Our envelopes feature a new tear-away thumb-notch opening that makes it easy to remove plates for processing. Made of soft, supple vinyl with rounded corners to optimize patient comfort, these envelopes are available in sizes 0, 1, 2 and 3.

HYGIENIC SLEEVES



2.5 x 4 cm



8.5 x 4 cm

Hygienic, single-use plastic sleeves to fit a variety of holders on the market. Akzenta Hygienic Sleeves are made of soft but strong plastic. They are elastic, latex-free, odourless and have no sharp seams. Tight fit with easy application and removal. Comes as a continuous roll of 500 sleeves within a dispensing box. Compatible with: Vatech Pax-i Digital OPG, Vatech Pax Primo Basic Digital OPG, Vatech Pax-Flex P Digital OPG, Vatech Pax-Flex 3D Panoramic X-Ray, Vatech Pax-Duo 3D CBCT & Panoramic X-Ray, Novus Digital OPG Machine, X-Mind Pano D+ Digital OPG Machine, Satelec X-Mind 3D Digital X-Ray. Bite Block Hygiene Sleeve 4.6 x 2.1cm Will fit the following OPG Machines: Vatech – Pax-Primo Basic Pax-l Pax-Flex P Pax-Flex 3D Pax-Duo 3D. Satelec – Novus X-Mind Pano D+ X-Mind 3D.



STERILISATION

TOP SAFE STERILISATION REELS

Disposable reels made of special paper and foil. Extremely safe due to 3 parallel running weld seams (multiline method).
Manufactured in accordance with EN 868. Strictly controlled production process.
Appropriated for steam sterilisation from 121° - 134° C, for gas sterilisation 30° – 54° C.
Successful sterilisation is indicated by the color changes of the indicator stripes. CE



Article Number	Format	Side fold	Rolls Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight Kg
SRAZ055	55 mm x 200 m	No	3	21 x 21 x 21	3.6	150	15	10	80 x 120 x 225	560
SRAZ075	75 mm x 200 m	No	2	21 x 21 x 21	3.2	150	15	10	80 x 120 x 225	500
SRAZ100	100 mm x 200 m	No	2	21 x 21 x 21	4.5	150	15	10	80 x 120 x 225	875
SRAZ150	150 mm x 200 m	No	1	21 x 21 x 21	3.3	150	15	10	80 x 120 x 225	695
SRAZ200	200 mm x 200 m	No	1	21 x 21 x 21	2.2	150	15	10	80 x 120 x 225	350
SRAZ250	250 mm x 200 m	No	1	21 x 21 x 41	3	75	15	10	80 x 120 x 220	245
SRAZ300	300 mm x 200 m	No	1	21 x 21 x 41	6.5	75	15	10	80 x 120 x 220	507

TOP SAFE STERILISATION POUCHES

Article Number	Size mm	Pieces Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Pallet Size CM	Pallet Weight KG
APSG-2-140260	140 x 260	6 x 200	31 x 31 x 20	6.92	81	80 x 120 x 195	580
APSG-2-190330	190 x 330	6 x 200	40 x 39 x 20	11.01	42	80 x 120 x 155	482
APSG-2-90250	90 x 250	9 x 200	31 x 31 x 22	6.36	81	80 x 120 x 213	535

Disposable sterilisation pouches. First-class paper- foil combination. For packing of the sterilisation goods. Autosealed sterilisation pouches with indicator.
Safe gas- and steam sterilisation based on effortless infiltration of steam and gas through the closed pouches.




TOP SAFE STERILISATION SHEETS

Article Number	Format cm	Color	Pieces Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
HY-52	40x40	White	500	65x63x7,5	5.9	30	-	-	80x120x130	197
HY-55	50x50	White	500	65x63x7,5	8.5	30	-	-	80x120x130	275
HY-56	60x60	White	500	65x63x7,5	11.03	30	-	-	80x120x130	351
HY-58	75x75	White	250	78x50x9	9.95	20	-	-	80x120x150	219
HY-61	80x80	White	250	93x50x9	11.1	20	-	-	80x120x150	242
HY-62	90x90	White	250	93x50x9	13.65	20	-	-	80x120x150	273
HY-63	100x100	White	250	105x55x8,5	16.6	20	-	-	80x120x152	352
HY-64	120x120	White	125	126x51x7	12.75	20	-	-	80x120x150	275
HY-67	40x40	Green	500	65x63x7,5	5.9	30	-	-	80x120x130	197
HY-71	50x50	Green	500	65x63x7,5	8.5	30	-	-	80x120x130	275
HY-72	60x60	Green	500	65x63x7,5	11.03	30	-	-	80x120x130	351
HY-73	75x75	Green	250	78x50x9	9.95	20	-	-	80x120x150	219
HY-75	80x80	Green	250	93x50x9	11.1	20	-	-	80x120x150	242
HY-76	90x90	Green	250	93x50x9	13.65	20	-	-	80x120x150	293
HY-77	100x100	Green	250	105x55x8,5	16.6	20	-	-	80x120x152	352
HY-78	120x120	Green	125	126x51x7	12.75	20	-	-	80x120x150	275
HY-79	140x140	Green	100	143x55x8	13.7	20	-	-	80x120x150	294

Article Number	Format cm	Color	Pieces Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
HY-108	0.5x100	White	2	15x15x10,1	6.93	30	-	-	80x120x110	228
HY-109	0.75x100	White	1	15x15x7,6	5.06	50	-	-	80x120x125	273
HY-110	1x100	White	1	15x15x101	6.93	30	-	-	80x120x110	228
HY-143	0.5x100	Green	2	23x34x51	6.93	30	-	-	80x120x110	228
HY-145	0.6x100	Green	1	15x15x61	4.16	50	-	-	80x120x125	228
HY-147	0.75x100	Green	1	15x15x76	5.06	50	-	-	80x120x125	273
HY-154	1x100	Green	1	15x15x102	6.93	50	-	-	80x120x110	228

60 g/qm
strict control of production
absolute hygienie
top quality
This sterilization paper features a water repelling texture.
Sheets produced in accordance with DIN 868-1 and DIN 868-2 chapter 2.14 and 2.15 for wrapping of sterilization goods.
Made in Switzerland.



TISSUE NAPKINS

EUZELLA TISSUE NAPKINS


Smooth, soft and absorbent.
Produced environmently friendly.
Hygienic packaging.
Color: White.



Article Number	Format cm	Number of layers	Pieces Pack	Packs Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
TISS-3-E214	40 x 40	2	250	6	57.5 x 40.5 x 20	7.44	32	4	8	80 x 120 x 175	258
TISS-3-E414	40 x 40	4	200	4	57.5 x 40.5 x 20	7.94	32	4	8	80 x 120 x 175	264

TOP SOFT Tissue Napkins

Smooth, soft and absorbent.
Produced environmently friendly.
Hygienic packaging.



Article Number	Format cm	Color	Number of layers	Pieces Pack	Packs Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
TISS-3-0632	40 x 40	White	2	250	6	57.5 x 40.5 x 20	7.44	4	8	32	80 x 120 x 175	258
TISS-3-0433	40 x 40	White	3	250	4	57.5 x 40.5 x 20	7.94	4	8	32	80 x 120 x 175	274
TISS-3-0434	40 x 40	White	4	200	4	57.5 x 40.5 x 20	7.94	4	8	32	80 x 120 x 175	274
TISS-3-0998	33 x 33	White	3	250	4	34 x 34 x 25	5.06	6	7	42	80 x 120 x 190	233
TISS-1-0384	33 x 33	Yellow	3	250	4	34 x 34 x 25	5.06	6	7	42	80 x 120 x 190	233
TISS-3-0406	33 x 33	Fresh Green	3	250	4	34 x 34 x 25	5.06	6	7	42	80 x 120 x 190	233
TISS-1-0203	33 x 33	Navy Blue	3	250	4	34 x 34 x 25	5.06	6	7	42	80 x 120 x190	264

TOOTHBRUSHES

TOP BRUSH



Yellow
ZABU2-AZ-Y



Orange
ZABU2-AZ-O



Pink
ZABU2-AZ-P



Lilac
ZABU2-AZ-L



Fresh Green
ZABU2-AZ-FG



Med Blue
ZABU2-AZ-MB

The TOP BRUSH has an integrated toothpaste with extra fresh taste.
Handy solid design, made of food safe and authentic polypropylene with transparent toothpaste



Dark Blue
ZABU2-AZ-B



Burgundy
ZABU2-AZ-WR



Fuchsia
ZABU2-AZ-F



Yellow, Orange,
Dark Green, Dark
Blue
ZABU2-AZ-MIX



Burgundy, Med Blue,
Fresh Green, Pink
ZABU2-AZ-MIX1



Yellow, Orange, Pink,
Pastel Blue
ZABU3-AZ-MIX

The TOP BRUSH PLUS toothbrushes with extra fresh taste and ultra-soft bristles made of nylon and oval brush head with transparent tooth paste. CE

Pieces Box	Boxes Carton	Carton Size CM	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
100	10	57 x 36.5 x 19.5	7.5	40	4	10	80 x 120 x 210	320

TRAY PAPER INSERTS

TOP SAFE PATIENT TOWELS 2-LAYER



White
TRAY18-N-W
TRAY28-N-W



Yellow
TRAY18-N-Y
TRAY28-N-Y



Orange
TRAY18-N-O
TRAY28-N-O



Pink
TRAY18-N-P
TRAY28-N-P



Burgundy
TRAY18-N-WR
TRAY28-N-WR



Lilac
TRAY18-N-L
TRAY28-N-L



Fresh Green
TRAY18-N-FG
TRAY28-N-FG

Highly absorbent and great strength against water,
ecologically friendly cellulose | 55gm²
Size: 18 x 28 cm and 28 x 36 cm



Blue
TRAY18-N-B
TRAY28-N-B



Green
TRAY18-N-G
TRAY28-N-G



Med Blue
TRAY18-N-MB
TRAY28-N-MB



Fuchsia
TRAY18-N-F
TRAY28-N-F



Size cm	Sheet Unit	Units Carton	Carton Size cm	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
18 x 28	250	10	37.5 x 29.5 x 23	7.3	64	8	8	80 x 120 x 198	487
28 x 36	250	5	37.5 x 29.5 x 23	7.3	64	8	8	80 x 120 x 198	487

WET CREPE NAPKINS

PREMIER | WHITE WET CREPE NAPKINS



White
NSKS-3-0110L

Strong and soft crepe paper.
Protects patients against moisture
and contamination.
Highly absorbent and great strength
against water, tear proof.
Size: 34 x 38 cm

TOP COVER | WHITE WET CREPE NAPKINS

White
NSKS4-AZ



Article Number	Papers Box	Boxes Carton	Carton Size cm	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
NSKS-3-0110L	1000	6	80 x 40 x 37	21.75	12	3	4	80 x 120 x 163	261
NSKS4-AZ	1000	4	72.5 x 39 x 17.5	13.8	27	3	9	80 x 120 x 172	392

Top Selection+

PREMIER PLUS | COLOURED + WHITE WET CREPE NAPKINS



White
CNH3610-N-W
CNH365-N-W



Yellow
CNH3610-N-Y
CNH365-N-Y



Green
CNH3610-N-G
CNH365-N-G



Blue
CNH3610-N-B
CNH365-N-B

Strong and soft crepe paper.
Protects patients against moisture
and contamination.
Highly absorbent and great strength
against water, tear proof.
Size: 36 x 36 cm
3610: 1000 pieces
365: 500 pieces

Papers Carton	Carton Size cm	Carton Weight Kg	Cartons Pallet	Cartons Layer	Layers Pallet	Pallet Size CM	Pallet Weight KG
1000	37 x 37 x 15	5.83	42	6	7	80 x 120 x 120	265
500	37 x 37 x 8	2.8	84	6	14	80 x 120 x 127	255

PRODUCT SETS

COCO COLOR SETS



Lilac
COCO-AZ-L



Fresh Green
COCO-AZ-FG



Burgundy
COCO-AZ-WR



Fuchsia
COCO-AZ-F



Pink
COCO-AZ-P



Orange
COCO-AZ-O



Med Blue
COCO-AZ-MB

Free Gift inserts in every carton!

1 package Dental Rolls of 300g
1 package Instant Cold Packs of 1 piece
1 package Micro Applicators of 4 cans

Carton Content:

Patient Towels of 50 pieces
Aprons of 80 pieces
Saliva Ejectors of 100 pieces
Cups of 50 pieces
Tray Papers of 250 pieces
Headrest Covers of 50 pieces
Face Masks of 50 pieces
Nitrile Gloves of 100 pieces in XS or S or M
Latex Gloves of 100 pieces in XS or S or M
Aspirator Tips of 10 pieces
Toothbrushes of 100 pieces

OUR PRODUCT LINES IN DETAIL

SMILY SET

We designed the Smily Sets to make sure our youngest patients enjoy every moment of their treatment with our fun and joyful Smily motives. We added a FREE special gift (Brushes or Stickers).

The Smily Set contains:

100pcs Smily Cups
100pcs Smily Masks
100pcs Smily Patient Towels

+ Free Gift:

100pcs Top Brush Plus (SMILY-SET1)
100pcs Emoji Stickers (SMILY-SET2)





TOOTHBRUSHES
Top Brush Plus

see page 55



FACE MASK
Top Mask Plus
Coolbreath

see page 24



APPLICATORS
Top Sticks Plus

see page 10

Designed With Excellence In Mind
The Top Selection+ delivers the most sophisticated selection of disposables. Aimed at those who demand nothing but the best, you will experience excellence like never before.



GLOVES
Top Touch Plus
Nitrile Ultrasoft

see page 28



BIB HOLDER
Top Easy Plus

see page 14



PROTECTIVE SLEEVES
Barrier Envelopes/
Hygienic Sleeves

see page 50

Fulfilling Your Everyday Needs

The Top Selection is the high quality line aimed at giving you superb disposables at very competitive prices. Experience quality and comfort in a new dimension.

GLOVES
Top Touch
Latex Powdered

see page 27



FACE MASK
Top Mask IIR
3 - Ply Non Latex

see page 23



CUPS
Top Cups

see page 17



HEADREST COVERS
Top Cover

see page 38



SALIVA EJECTORS
Top Flex

see page 41



BIB HOLDER
Top Easy

see page 14



PATIENT TOWELS
Top Safe

see page 44



LET COLOR TOUCH YOU

*Make Your Dental
Practice Colourful
and Your Patients
Smile*



OUR 11 MOST POPULAR PRODUCTS IN 7 VIBRANT COLORS.

Strengthen your practice identity with visual impact with your favorite color.

The COCO COLOR COLLECTION is a carefully selected line of products that are in constant demand by dentist and medical offices across the world.

We make you shine

With the COCO COLOR COLLECTION, Akzenta is giving dental and medical practices a fresh and consistent look and feel. All 11 products of our COCO COLOR COLLECTION can be bought in one of 7 vibrant colors that are matching perfectly throughout the whole product line.

VIBRANT COLORS

Choose from one of 7 colors and equip your dental or medical practice with our complete set of color-matching, premium products.

-  MED BLUE
-  LILAC
-  PINK
-  FUCHSIA
-  LAGUNA BLUE
-  ORANGE
-  FRESH GREEN

see page 58

SMILY Series

Making your little patients smile

The Smily Series aims at making visits to the doctor's office more joyful for children. Our cute design helps putting children at ease and create a friendly atmosphere.

SMILY PATIENT TOWELS

The Smily Patient Towels are highly absorbent (on top) and waterproof (back side) to ensure a high comfort of your patients.
see page 48



SMILY MASK

The Smily Mask is distinguished by its high wearing comfort even during extended wearing time.
see page 25



SMILY CUPS

The high quality Smily Cup is made of food safe, authentic polypropylene: comfortably smooth, yet very resistant and solid.
see page 21



SMILY SETS

Akzenta introduces It's new Smily Sets, building upon the success of the very popular Smily Series.
see page 58

BLACK Series



TRY SOMETHING NEW AND EXCITING!

The Black Series is aimed at professionals that want to showcase their personality by using our pitch black products. Make a statement - with the Black Series by Akzenta.



see page 32



see page 25

LATEX BLACK

NITRILE BLACK



Swiss Edition



Swiss Excellence For Your Satisfaction



FACEMASKS

see page 22

PATIENT TOWELS

see page 44

Pure Edition



Our bestselling Pure Editon Gloves are the softest gloves you will find. Its silky smooth surface makes long wearing periods a breeze.



LATEX SENSITIVE

see page 33

NITRILE SOFT

see page 33

Paradise Edition



Introducing the attractive motifs of our paradise edition



PARADISE CUPS

see page 20



PARADISE TOWELS

see page 25



PARADISE MASKS

see page 25

GET TO KNOW US

GET IN TOUCH



SALES

Department
tel. +41 91 210 46 08
fax. +41 91 682 85 20

e-mail: sales@akzenta.com



CUSTOMER SERVICE

Department
tel. +41 91 210 46 06
fax. +41 91 682 85 20

e-mail: service@akzenta.com



AKZENTA INTERNATIONAL SA

Via Giuseppe Motta, 24
CH-6830 Chiasso

tel. +41 (0) 91 921 14 92
fax +41 (0) 91 682 85 20

e-mail: business@akzenta.com

www.akzenta.com



Scan for vCard

PURCHASING

Department
tel. +41 91 210 46 05
fax. +41 91 682 85 20

e-mail: purchasing@akzenta.com



MARKETING

Department
tel. +41 91 210 46 09
fax. +41 91 682 85 20

e-mail: marketing@akzenta.com



ACCOUNTING

Department
tel. +41 91 210 46 03
fax. +41 91 682 85 20

e-mail: accounting@akzenta.com



WE ARE LOOKING FORWARD TO HEARING FROM YOU!

We are a motivated and innovative team located in the beautiful Tessin in the Italian part of Switzerland. Due to our efficient company structure we are able to react to changing markets and customer demands fast and reliably. We are always open for your ideas and wishes. Having a personal conversation with you is of utmost importance to us. Constantly on the move we are always at your service locally. We treasure your contact!



AKZENTA INTERNATIONAL SA

Via Giuseppe Motta, 24
CH-6830 Chiasso

tel. +41 (0) 91 921 14 92

fax +41 (0) 91 682 85 20

e-mail: business@akzenta.com

www.akzenta.com

